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Pedalo Federbrett Vestimed 50 inkl. Textil-Powerband

Art.-Nr. 02309

▶▶▶ [zum Produkt ...](#)

Pedalo Trainer

▶▶▶ [zur Kategorie...](#)

Anleitung

14004050

Pedalo® - Vestimed 50

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English

Patent 102007048399. The special feature of **Vestimed** is the spring system, which allows all exercises to be carried out in an instable stand. The impulses to the nervous system and muscles are intensified and the training becomes more effective. Posture and movement coordination will be improved noticeably. **Vestimed®** is a total body training device. It is suitable for all ages and enhances individual skills of your brain, nervous system, arms, legs, back, trunk, stomach. The training has a targeted effect on improving the coordination skills. Training is done while standing up, sitting or in quadruped position and with or without the use of hand and foot straps. Hand and foot straps are quickly mounted per carbine. **Vestimed®** can be used either for leisure, health, or pro-sports, as well as in therapy after injuries to the loco-motor system and for people with neurological deficiencies as in Parkinson. Preselected spring strengths, adapted to body weight, provide optimal effect. By the use of the new hooking-in system, an individual pre-tension of hand and leg straps, adapted to users' physical capability, can be adjusted.

Item no. 140 04050

Specification:

Senso Springs with optimal effect from 30-110 kg load, platform ø 50 cm, 18 cm high, birch plywood, 1 pair of Textile-Power Bands washable with loops, 2 fixing loops, 6.7 kg.

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English

Body perception

a) Stand with both feet on the Vestimed® and notice the movements of your body. Quit your tension completely and play with your balance.

b) Seesaw forth and back gently and relaxed by shifting weight between the forefoot and heel. Then seesaw to the left and right by weight shift to the left leg and to the right leg.

c) Align your body in an upright posture and hold the position.

d) Switch to tiptoe-position and hold the position

e) Go to slight crouch through to downhill position and hold the position.

If you feel safe, perform the movements with your eyes closed.

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English

Visual perception

a) Stand with both feet on the Vestimed® and notice the movements of your body. Look at a fixed point at eye level directly in front of you.

b) Based on a) turn your head to your left shoulder, find your balance, then turn your head to your right shoulder, find your balance, then turn back to the middle and find your balance again. Alternatively with outstretched arms and a ball in one hand.

c) Based on a) turn your head in short quick movements several times back and forth (like saying “no”). Find your balance.

d) Based on a) move your head in circle, find your balance.

If you feel safe, perform the movements with your eyes closed.

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English

Motion training


a) Stand with both feet on the Vestimed® and notice the movements of your body. Quit your tension completely and play with your balance.

b) Step into a regular rhythm on the spot. The more secure you feel, the longer you stay on one leg.

c) As b) in addition swing your arms.

d) As c) change also your direction of view.

If you feel safe, perform the movements with your eyes closed.

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English

Stair climbing

a) Get alternately slowly and equally with one leg on the Vestimed® and each back down again. Gradually increase the time spent on one leg on the Vestimed®.

b) As a), proceed in a series of steps over the Vestimed®. While descent pay attention to do it slow.

If you feel safe, perform the movements with your eyes closed.

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English

Leg axis stability

a) Stand on one foot centered on the Vestimed® and keep your balance.

b) Play with your free leg. Stretch it alternating times forward, sideways, backwards.

c) As a) move your foot in addition in a circle, change your direction of view.

If you feel safe, perform the movements with your eyes closed.

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English

a) Lateral trunk stabilization

- 1) Go to the lateral forearm support and stabilize a straight posture.
- 2) Stretch your free leg and arm in the extension of the body, then run elbow and knee to each other.

b) Push-up position

- 1) Go into push-up position and keep only the upper position.
- 2) Go into push-up position and keep only the lower position.
- 3) Perform several active pushups.

c) Forearm support

- 1) Go into forearm support with your feet on the Vestimed®.
 - a. Stabilize.
 - b. Alternately lift one leg slightly.
 - c. Stretch alternately one hand forward.
- 2) As 1) but forearm support on Vestimed®.

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English

Diagonal pelvic stabilization:

a) Stand in stepping position (one leg is on the Vestimed®). Perform your rearward knee toward the floor and stabilize your pelvis while doing this.

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Exercises with the Vestimed® + Powertubes

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English

Core training

- a)** Go into a slightly squat and hold the power bands on tension in front of your body. Perform squats.
- b)** As a) but lead your crossed arms at the same time in front of your body from the bottom up against the resistance of the power bands until they have reached their full stretching.
- c)** Lead your arms alternately in diagonal against the resistance from bottom to top.
- d)** as b) but but when you drag with the right arm, the left leg is also splayed (one leg stance) and vice versa.

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English

Training of shoulder girdle

a) Stabilize your arms and shoulder girdle against the resistance of the power bands in different posture positions.

- Arm is stretched forward / sideways

b) As a) but one arm is stabilized and the second arm is led in different positions from bottom to top.

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English

Raise of the upper part of your body

- a)** Stand with both feet on the Vestimed® (alternatively one-legged) and hold the power bands in front of you.
- b)** As a) but carry the power bands from the front to the top in stretching.
- c)** As a) but one arm is still in front of your body (alternately move one arm)
- d)** As you pull, stretch alternately one leg backwards. While doing this, stabilize your pelvic.

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English

Full balance training

One foot is secured in the attachment loop so that you can work with resistance against the power band. Depending to the exercise on the Vestimed®, place the device as the power band is mounted in line for movement execution.

a) Move in the one leg stand your free leg against the resistance forward, backward, sideways from inside to outside (abduction).

b) As a) but you keep with the free leg opposite hand the power band resistance.

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
English

a) Shoulder girdle training

Stretch from the supine position the arms, which are first outstretched towards the pelvis, against the resistance of the power bands over your head.

b) Stabilization of back stretching

as a) additionally lift the pelvis that knee - pelvis - shoulder form a straight line.

c) Abdomen-trunk-training

Hold power bands shortly, lift your upper body towards your knees, hold position for a short time, Slowly quit the tension.

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Training concepts with the Vestimed®

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English

naturalMOVE® are training concepts for everybody, which connects natural and psychological forms of motion, exercises for muscular functions and Core-Training in a unique way. The heart of the **naturalMOVE®** systems are the 3D balance training systems of pedalo®. The special training concept offers several fields of application and reaches both „the young audience“ and seniors. The training successes are a result from the unique combination of functional total body workout with coordinative demands. Through specialist trainers seminars, which are patented by pedalo® Academy, physical therapists and health fitness trainers are with a single training able to sustainably implement several training concepts.

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Dear customer,

The pedalo®-article you have chosen is made of natural wood. Differences in colour, structure, grain, small included knobs or scars prove the genuineness of our product. They are natural and therefore features of using original wood and no quality reduction.

Please note:

For the use of this device a certain degree in motor abilities and circumspection is assumed. In case of doubt first use device only with assistance.

copyright by Holz-Hoerz GmbH, printed copies also partially are not allowed.

==== www.sport-tec.de ====

Sport-Tec GmbH
Physio & Fitness

Lemberger Straße 255
D-66955 Pirmasens

Tel.: +49 (0) 6331 1480-0

Fax: +49 (0) 6331 1480-220

E-Mail: info@sport-tec.de

Web: www.sport-tec.de