DE 

EN 

FR 

IT 

ES 

NL 

SPORT-TEC ARM- UND BEINTRAINER MOVE 5.0

SPORT-TEC MINI EXERCISE BIKE MOVE 5.0

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Ihr Sport-Tec Team

WELCOME TO SPORT-TEC

Thank you for your purchase. We are pleased that you have opted for a Sport-Tec product and hope that you really enjoy using it. In this manual, you will find all the important information for use. If you have any further questions, please contact us – we will be happy to help you.

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S. 4-9

EN - MANUAL

Read this manual carefully before use. Follow the instructions, and keep this user manual in a safe place.

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
GENERAL INFORMATION

The Mini Exercise Bike move 5.0 is perfect for your home workout. Make yourself comfortable and train your arms and legs with just one training device. You can vary the resistance using the adjusting wheel. Thanks to the simple operation, this training device is ideal for older people. The drive belt allows fluid movements and precise regulation of the resistance.

SCOPE OF DELIVERY

No.	Designation	Quantity
1	Main Frame	1x
2	Front Stand	1x
3	Rear Stand	1x
4	Pair of pedals (right/left)	1x
5	Washer \varnothing 6 mm	4x
6	Screw M6x15	4x
7	AAA Battery	1x
8	Wrench	1x
9	Allen Key	1x

SAFETY INSTRUCTIONS

To make sure that you do not have any health or physical limitations that affect your exercise, consult your doctor before starting the exercise routine. This is particularly important when taking medication. In addition, always pay attention to the signs your body communicates and stop exercising immediately if symptoms such as chest pressure or irregular heartbeat occur, or if you feel pain or feel uncomfortable. In this case, contact your doctor immediately.

Before each use, check that the device is in good condition and make sure that the screws are tightened. In the event of defective components or unusual noises during use, interrupt the workout and correct the problem before continuing the training.

Make sure that the Mini Exercise Bike is placed on a stable and flat surface. Protect your floor with a mat if necessary.

When using the device, wear suitable clothing that offers sufficient freedom of movement and cannot get caught in the device.

The adjusting wheel for regulation of the resistance should be set to a low level at the beginning of the training. Check the resistance by turning the pedals. If necessary, consult your doctor or therapist regarding the choice of resistance.

Never stand up while using the Mini Exercise Bike. This can lead to injuries and damage the device. Exercise only while sitting.

ASSEMBLY

Important instructions:

Before starting the montage, first read through all the steps carefully.

Preparation:

Unpack all parts and lay them out clearly on the floor.

First, check the condition and completeness of all parts. Never use parts that have already been damaged. Contact us in the event of damage or incompleteness.

Stand montage

Attach the front stand **2** and the rear stand **3** to the main frame **1**. Use the washer **5** and screws **6** provided.

Pedals installation

Attach the left pedal **4** to the left crank and the right pedal **4** to the right crank. To attach, turn the right pedal clockwise and the left pedal counterclockwise.

NOTES ON USE

Display

Switch on the display by pressing the red button.

Resistance

The further you turn the adjustment wheel, the greater the resistance.

Adjust pedal loops

The pedal loops are adjustable in two sizes. They can be loosened and repositioned using a screwdriver.

DISPLAY

The display turns on as soon as the pedals are moved. You can manually switch on the display by pressing the red button.

After 4 minutes of inactivity (no movement is detected) the display switches off automatically.

The following information can be displayed:

- Training time (T),
- Distance (D),
- Completed revolutions (C),
- Total revolutions (T-C),
- Calories burned (CA),
- Training time, distance, completed revolutions, total revolutions and calories burned alternately (SCAN).

By pressing the red button in succession, you can activate / deactivate SCAN-mode and switch between the different display modes.

RESET DISPLAY

To reset all values, press and hold the red button for 4 seconds. The total revolutions (T-C) will only reset when the battery voltage is interrupted.

TRAINING TIME (T)

To display the training time, press the red button below the display until the arrow (▲) points to „T“ but „SCAN“ is not flashing. If „SCAN“ is flashing, all information is displayed alternately.

Display format: 00:00 - 99:59 (HR:MIN)

DISTANCE (D)

To display the distance you have completed, press the red button below the display until the arrow (▲) points to „D“ but „SCAN“ is not flashing. If „SCAN“ is flashing, all information is displayed alternately.

Display format: 0.000 - 0.999

COMPLETED REVOLUTIONS (C)

The mini exercise bike move 5.0 counts the revolutions. To display the completed revolutions, press the red button below the display until the arrow (▲) points to „C“ but „SCAN“ is not flashing. If „SCAN“ is flashing, all information is displayed alternately.

Display format: 0 - 9999

TOTAL REVOLUTIONS (T-C)

The mini exercise bike move 5.0 counts the revolutions. To display the total revolutions, press the red button below the display until the arrow (▲) points to „T-C“ but „SCAN“ is not flashing. If the „SCAN“ is flashing, all information is displayed alternately.

Display format: 0 - 9999

CALORIES BURNED (CA)

To display the calories burned, press the red button below the display until the arrow (▲) points to „CA“ but „SCAN“ is not flashing. If the „SCAN“ is flashing, all information is displayed alternately.

Display format: 0 - 9999.9

ALTERNATING DISPLAY (SCAN)

To display the training time, the distance, the completed revolutions, the total revolutions and the calories burned alternately, select a display option and make sure that „SCAN“ is flashing at the same time. The information is now displayed alternately. The arrow (▲) indicates which information is currently displayed.

CHANGE BATTERY

If you have problems with displaying the information clearly replace the battery of the display. To do this, remove the display by gently pulling on the frame. The battery is on the back. Make sure that you insert the battery correctly.

CARE AND CLEANING

Switch off the device, then clean using a non-abrasive cleaner or mild detergent with a soft cloth. Abrasive cleaners and / or abrasive cleaning pads may seriously damage the product beyond repair and should not be used.

Always ensure that you wipe dry the equipment after cleaning.

DISPOSAL

You must not dispose of this device in domestic waste. Most local authorities have a dedicated collection system for electrical equipment and disposal is free for end users.

PRODUCT FEATURES

	Sport-Tec Mini Exercise Bike move 5.0
Item no.	22004
Colour	black / white
Battery	1,5 V AAA
Weight	3,5 kg
Dimensions (LxBxH)	39x33x30 cm
Country of origin	Made in China

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Sport-Tec GmbH