DE 

EN 

FR 

IT 

ES 

NL 

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Ihr Sport-Tec Team

WELCOME TO SPORT-TEC

Thank you for your purchase. We are pleased that you have opted for a Sport-Tec product and hope that you really enjoy using it. In this manual, you will find all the important information for use. If you have any further questions, please contact us – we will be happy to help you.

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S. 4-11

EN - MANUAL

Read this manual carefully before use. Follow the instructions, and keep this user manual in a safe place.

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GENERAL


The move 8.0 is a motorised arm and leg trainer. It enables passive training and is intended to help strengthen the muscles.

The device is ideal for your home workouts. The short training sessions can be integrated into everyday life and are carried out in comfort while seated.

The removable remote control allows you to control the move 8.0 conveniently. You can choose between 12 speed levels, the directions forwards and backwards, and different time intervals. What's more, you have a total of 6 training modes available so you can carry out your specified training.

SCOPE OF DELIVERY

No.	Designation	Quantity
1	Main body with remote control	1x
2	Footboard	2x
3	Screw	4x
4	Right pedal	1x
5	Left pedal	1x
6	Spanner	1x

SAFETY INSTRUCTIONS

Please note that this device has been developed exclusively for home use.

Improper use of the arm and leg trainer can lead to injuries.

Before each use, check the proper condition of the device and make sure that the screws are securely tightened. If any components are defective or if there are any abnormal noises during use, stop exercising immediately. If necessary, have the device checked by qualified personnel before using it again.

Make sure that the mains supply voltage is suitable for the device before plugging it in.

Do not use the device if the mains cable becomes worn or damaged. Keep the mains cable away from hot surfaces and do not place it under a carpet or rug. Make sure that the mains cable is not in the way and that no one can trip over it.

Do not permit children to play with this device.

Never leave the device unattended when it is plugged in. Unplug it from the wall socket after each exercise session. To unplug the device first put the switch to the "0" position and then unplug from the wall socket.

It is not advisable to keep the device in a damp place for long periods as it will be prone to rust.

Do not use the device in bathrooms, or in wet or damp locations, or outdoors.

Please consult a medical professional before using the device for the first time or before undertaking a new form of exercise. Do not use the device if you feel unwell or immediately after eating. Always pay attention to your body's signals during training and stop training immediately if you experience symptoms such as pressure in your chest, an irregular heartbeat, pain or discomfort. In this case, contact your doctor immediately.

Make sure that you warm up sufficiently before training and plan time to cool down slowly after training.

Place the device on a flat surface, with at least 0.5 metres of free space at the sides and rear, for exercise. Wear suitable clothing and footwear for exercising that cannot get caught in the pedals. Ensure that the cable does not come into contact with the pedals while in use.

Do not touch or insert any objects into the device while the device is in use.

Keep the removable remote control handy during training so that you can intervene in the training at any time and, if necessary, stop it immediately.

You can stop training at any time by pressing the red  button on the remote control.

Never stand up on the arm and leg trainer. This can result in injuries and damage to the device. Exercise only while seated.

Warning! Risk of electric shock and injury!

Never carry out repairs on the device yourself. We recommend that repairs to electrical devices only be carried out by the manufacturer or the responsible customer service department.

DE

EN

FR

IT

ES

NL

BEFORE USE

Carefully remove all packaging. Avoid using knives or other sharp objects, as this may damage the surface of the product.

Check the device for any visible damage. If you see any damage or suspect a fault, do not use your device, but instead contact our service for support.

ASSEMBLY

Important notes:

Before assembling the move 8.0, first read through all the instructions carefully.

Adjust the left-hand pedal by turning it COUNTER clockwise. Adjust the right-hand pedal by turning it clockwise. Adjusting it the wrong way will damage the thread and invalidate your warranty.


Note:

The words right and left refer to the position that the user adopts when exercising.

PREPARATION


Unpack all parts and lay them out clearly on the floor.

First, check the condition and completeness of all parts. Never use parts that have already been damaged. Contact us if any parts are damaged or missing.

Mounting the feet

Attach the foot parts ② to the main body ① with the screws ③ (x4) provided.

Attaching the pedals

Insert the right-hand pedal ④ onto the right-hand crank. Note the R marking on the pedal and crank. Make sure that you turn the pedal clockwise.

Insert the left-hand pedal ⑤ onto the left-hand crank. Note the L marking on the pedal and crank. Make sure that you turn the pedal counterclockwise.


Note: Make sure that you insert the pedals slowly to avoid incorrect screwing and damage to the thread.

MAIN SWITCH

To turn on the device, plug in the power plug and then move the main switch to the "1" position.

To switch off the appliance, first set the main switch to the "0" position and then disconnect the mains plug from the socket.

Never reverse the procedure, as the device may be damaged.

REMOTE CONTROL WITH DISPLAY

Press the POWER button on the remote to turn the device on.

If the arm and leg trainer is not operated within 10 seconds after turning on, the screen will go into sleep mode.

SPEED

Press the speed +/- buttons to increase or decrease the speed. Select speeds levels 1 to 12.

TIME

Press the time button +/- to select the desired training time: 3, 6, 9, 12 or 15 minutes.

When choosing the exercise mode, the time that appears by default is 3 minutes. If you want to increase your exercise time, you must do so after selecting the program, but before starting the exercise. The time cannot be modified once the exercise has started.

POWER

Training can be paused by briefly pressing the POWER button. When you are ready to resume the workout, briefly press the POWER button again.

MODE

When turning on the arm and leg trainer, the mode "0" is set – this is manual mode. You can change speed, time and direction in this mode. To choose from the 6 available auto programs, keep pressing the mode button till you reach your desired program.


DIRECTION

Press the direction button to switch between forward and backward motion.

STOP

When you press the STOP button, the arm and leg trainer turns off automatically.

DE

EN

FR

IT



ES

NL

ACCESSORIES

NON-SLIP MAT


To prevent the arm- and leg trainer from slipping, we recommend placing the non-slip mat underneath. This enables you to train safely and comfortably.

TRAINING


POSITION 1: SEATED (LOWER BODY)

- Use a stable, sturdy chair.
- Place the device at a comfortable distance from the chair to avoid overstretching your legs.
- Make sure that your lower back is supported. We recommend that you place a small cushion at the base of your back.
- Sit on the chair and place your feet on the pedals.
- Remove the remote control from the housing of the device. Turn on the device and select the desired training program to start your workout.

POSITION 2: SEATED (UPPER BODY)

- Use a raised, solid surface such as a table so that the device is at chest height with the side with the remote control facing you. Sit on a stable chair in front of the table.
- Place the unit at a comfortable distance to avoid overstretching your arms.
- Turn on the device and select the desired training program to start your workout.

DESCRIPTION OF MODES

TIME INTERVALS

Total training duration (min.)	3	6	9	12	15
Time per interval	18 sec.	36 sec.	54 sec.	1 min. 12 sec.	1 min. 30 sec.

THE MODES

MODE 1	Speed level	1	1	1	1	1	1	1	1	1	1
	Raotation direction	MANUAL									

Enter the level of speed (SPEED) and time (TIME), using the buttoms (+ and -) on the remote control. Once the training has started, you can change the speed (+/- buttons) and the pedaling direction (DIRECTION button) during the exercise, but you cannot change the set time. The time can only be modified BEFORE starting the exercise.

MODE 2	Speed level	1	2	3	4	4	4	4	3	2	1
	Raotation direction	MANUAL									

Enter the desired training time (TIME), using the buttons (+ and -) on the remote control. Once the training has started, you will not be able to modify the selected time, but you can change the speed (+/- buttons) and the pedaling direction (DIRECTION button) during the exercise.

If you modify the speed, the following sections will increase proportionally. For example, if during the first section you increase the speed from level 1 to level 4 (+3 levels), when it goes to the second section the level will be 5, in the third section 6, and so on, adding 3 levels to what is indicated in the table above in each of the sections.

MODE 3	Speed level	1	4	1	4	1	4	1	4	1	4
	Raotation direction	MANUAL									

Enter the desired training time (TIME), using the buttons (+ and -) on the remote control. Once the training has started, you will not be able to modify the selected time, but you can change the speed (+/- buttons) and the pedaling direction (DIRECTION button) during the exercise.

If you modify the speed, the following sections will increase proportionally. For example, if during the first section you increase the speed from level 1 to level 4 (+3 levels), when it goes to the second section the level will be 7, in the third section 4, and so on, adding 3 levels to what is indicated in the table above in each of the sections.

MODE 4	Speed level	1	1	1	1	1	1	1	1	1	1
	Raotation direction										

Enter the desired training time (TIME), using the buttons (+ and -) on the remote control. Once the training has started, you will not be able to modify the selected time, but you can change the speed (+/- buttons) during the exercise. In this program, the speed is constant during all the sections, so if you modify it at some point during the training, it will remain at the new selected level until the end of the time.

For example, if during the first section you increase the speed from level 1 to level 4 (+3 levels), the speed will remain at level 4 during all the sections until the end of the selected time. In this program you will not be able to modify the pedaling direction either, which will automatically change from the pre-established way, alternating in each of the 10 equal time intervals into which the entire selected time is divided. For example, if a time of 3 minutes has been selected, we will have 10 intervals of 18 seconds each. In another example, if the selected time were 15 minutes, we will have 10 intervals of one minute and a half each of them.

MODE 5	Speed level	1	1	4	1	4	1	4	1	4	1
	Raotation direction										

Enter the desired training time (TIME), using the buttons (+ and -) on the remote control. Once the training has started, you will not be able to modify the selected time, nor the pedaling direction, which is already preconfigured, but you can change the speed (+/- buttons).

If you modify the speed, the following sections will increase proportionally. For example, if during the first section you increase the speed from level 1 to level 4 (+3 levels), when it goes to the second section the level will be 4 still, in the third section it increases to 7, and so on, adding 3 levels to what is indicated in the table above in each of the sections. The pedaling direction will automatically change from the pre-set way, alternating in each of the 10 equal time intervals into which the entire selected time is divided. For example, if a time of 3 minutes has been selected, we will have 10 intervals of 18 seconds each. In another example, if the selected time were 15 minutes, we will have 10 intervals of one minute and a half each of them.

MODE 6	Speed level	1	1	1	1	1	1	1	1	1	1
	Raotation direction	CHANGES EVERY 5 SECONDS									

Enter the desired training time (TIME), using the buttons (+ and -) on the remote control. Once the training has started, this program will not allow you to make any changes to any of the parameters time, speed or direction, all of them are pre-established. The speed in this exercise mode will always be constant, at level 1, and the direction of rotation will alternate, changing every 5 seconds for the entire duration of the program.

CLEANING

Power off and unplug the device. Then clean the unit with a nonabrasive, mild detergent using a soft cloth. Abrasive cleaners and/or abrasive cleaning pads may seriously damage the product beyond repair and should not be used. Always ensure that you wipe dry the device after cleaning and that no water/liquid enters the interior of the device.

CARE AND MAINTENANCE

Please carry out a safety check of the product at regular intervals or if you have any concerns. Pay attention to fastening elements and tighten if necessary. If in doubt about the safety of the device, do not use it. Please contact our service if you require technical assistance.

DISPOSAL

Do not dispose of this device with household waste. Most local authorities have a special collection system for electrical appliances. Disposal is free of charge for the end user.

PRODUCT FEATURES

	Sport-Tec Motorized Mini Exercise Bike move 8.0 for arms and legs, with display
Art.-Nr.	22006
Colour	05-black-grey / 06-white-black
Power supply	Mains plug
Nominal voltage	AC 220-240 V
Rated frequency	50/60 Hz
Performance	70 W max.
Class	HC
Max. user weight	120 kg
Weight	4,5 kg
Dimensions (LxBxH)	42x40x29 cm
Country of origin	Made in China

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Sport-Tec GmbH