DE 

EN 

FR 

IT 

ES 

NL 

## SPORT-TEC ARM- UND BEINTRAINER MOVE 1.0

SPORT-TEC MINI EXERCISE BIKE MOVE 1.0

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---

## WELCOME TO SPORT-TEC

Thank you for your purchase. We are pleased that you have opted for a Sport-Tec product and hope that you really enjoy using it. In this manual, you will find all the important information for use. If you have any further questions, please contact us – we will be happy to help you.

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---

## EN - MANUAL

Read this manual carefully before use. Follow the instructions, and keep this user manual in a safe place.

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## GENERAL INFORMATION

The Mini Exercise Bike move 1.0 enables training at home or in the office. The Mini Exercise Bike is designed for light movement training. The resistance is infinitely adjustable.

## SAFETY INSTRUCTIONS

Please note that this device has been developed exclusively for home use.

Improper use of the arm and leg trainer can lead to injuries.

To make sure that you do not have any health or physical limitations that affect your exercise, consult your doctor before starting the exercise routine. This is particularly important when taking medication. In addition, always pay attention to the signs your body communicates and stop exercising immediately if symptoms such as chest pressure or irregular heartbeat occur, or if you feel pain or feel uncomfortable. In this case, contact your doctor immediately.

Before each use, check that the device is in good condition and make sure that the screws are tightened. In the event of defective components or unusual noises during use, interrupt the workout and correct the problem before continuing the training.

Make sure that the Mini Exercise Bike is placed on a stable and flat surface. Protect your floor with a mat if necessary.

When using the device, wear suitable clothing that offers sufficient freedom of movement and cannot get caught in the device.

The adjusting wheel for regulation of the resistance should be set to a low level at the beginning of the training. Check the resistance by turning the pedals. If necessary, consult your doctor or therapist regarding the choice of resistance.

**Never stand up while using the Mini Exercise Bike. This can lead to injuries and damage the device. Exercise only while sitting.**

## PRODUCT FEATURES

	Sport-Tec Mini Exercise Bike move 1.0
Item no.	22331
Colour	black / silver
Weight	2 kg
Dimensions (LxBxH)	50x42x28 cm
Country of origin	Made in China

## NOTES ON USE

1. Set the adjustment wheel to the desired resistance. Turn the adjustment wheel clockwise to increase the resistance and counterclockwise to reduce the resistance.

**i** You can pedal forwards or backwards to work different muscle groups.

2. The pedal loops are adjustable in two sizes. They can be loosened and repositioned using a screwdriver. For a better grip and to prevent slipping, the pedal loops should enclose your feet.

### APPLICATION ON THE FLOOR

1. Place the exercise bike in a comfortable position in front of you. If necessary, place a mat underneath to protect the floor.
2. Take a seat and put your feet on the pedals. Now start pedaling.

### APPLICATION ON THE TABLE

1. Place the exercise bike on the table top and sit on a chair in a comfortable position in front of it.
2. Keep your back straight. Grasp the pedals with your hands and start cranking. Try not to lean forward.

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**Sport-Tec GmbH**